


In Japan, Agency for Natural Resources and Energy is considering to amend the law “Enforcement Regulation of controlling the quality of volatile oils”, which is to strictly control the quality of biodiesel fuels (BDF).

As a result, a regulation of the compounds shown in the right table has been added. FAME (Fatty acid methyl esters) and TG (Triglycerides) in BDF are required to be analyzed using the HPLC method. Concentrations of FAME and TG are obtained from their calibration curves that have been already prepared using methyl stearate and trilinolein as their standard, respectively.

(K.Suzuki)

Target	Method
FAME (Fatty acid methyl esters)	HPLC-RI
TG (Triglycerides)	HPLC-RI
Methanol	GC-oxygen detection
Acid Number (Value)	Potentiometer Titration
Formic acid, acetic acid, propionic acid	Ion chromatography
Oxidation stability (Increased volume of acid number)	Potentiometer Titration

A chromatogram obtained from standard solution



Conditions

Column : Inertsil SIL-100A
(5 μm, 250 x 4.6 mm I.D.)
Column Cat. No. : 5020-01712
Eluents : A) n-Hexane
B) 2-Propanol
A/B = 996 / 4, v/v
Flow rate : 1.0 mL/min
Column Temp. : 40 °C
Detector : RI
Injection volume : 10 μL

Column performance test

In the HPLC method, following requirements should be fulfilled.


(Both values should be examined with 5- μL injection.)

- The retention time of methyl stearate is more than 3.5 minutes.
- Resolution factor (Rs) between the methyl stearate peak and trilinolein peak is more than 3.

Results

Retention time of methyl stearate: **4.05 min**

Resolution factor: **4.7**




Sample No.	Concentration (μg/10 μL)	
	Methyl stearate	Trilinolein
1	10	-
2	100	1
3	200	2
4	500	5
5	1000	10

Calibration curves and concentrations of standard solutions

A chromatogram obtained from BDF

The sample solution was injected after filtration using a 0.45- μ m membrane filter and dilution with light oil.




Visit our Website on

<http://www.glsciences.com/>

<http://www.inertsil.com/>

GL Sciences also offers an extensive range of GC and LC consumables.



Contact us

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-chome, Shinjuku-ku, Tokyo, 163-1130 Japan

TEL: +81-3 (5323)6620 FAX: +81-3 (5323)6621

GL Sciences, Inc. USA

4733 Torrance Blvd. Ste 255, Torrance, CA 90503

Tel: (310)265-4424 FAX (310)265-4425

Distributors Outside of Japan and USA

GL Sciences uses distributors in many countries.

You can find a local distributor in your country in the following url.

<http://www.glsciences.com/products/contact.html>

E-MAIL: world@glsc.co.jp

