


# Analysis of Putrefactive Non-Volatile Amines in Food by Pre-column HPLC

This note describes a determination method for non-volatile putrefactive amines in food using pre-column derivatization.

Putrefactive amines in food are produced by microorganisms. Many of the biogenic amines are known to be responsible for food-poisoning. To determine these amines, HPLC systems coupled with pre-column or post-column derivatization methods are frequently used.

Based on the inspection guideline by Japan Food Hygiene Association, pre-column derivatization method using DNS-Cl (5-dimethylaminonaphthalene-1-sulfonyl chloride) was adopted in this note. The putrefactive amines were successfully fluorescently derivatized, and the calibration curves showed excellent linearity. (K.Suzuki)


## Chemical structures



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovayion Software, Inc.

## Conditions

- Column Inertsil ODS-SP (5  $\mu$  m, 250 x 4.6 mm I.D.)  
Cat.No. 5020-02746
- Guard column Cartridge Guard Column E ODS-SP (5  $\mu$  m, 10 x 4.0 mm I.D.)  
Cat.No. 5020-08520
- Eluent A) CH<sub>3</sub>CN B) H<sub>2</sub>O  
A/B= 65/35, v/v, 1.0 mL/min
- Column Temp. 40 °C
- Detection FL Ex 325 nm, Em 525 nm
- Injection Vol. 10  $\mu$  L




The calibration curves for the putrefactive amines

## A chromatogram obtained from standard solution

Standard solution derivatized by DNS-Cl was analyzed.

1. Putrescine (5.0 mg/L)
2. Cadaverine (5.0 mg/L)
3. Histamine (100 mg/L)
- I.S. 1,8-Diaminooctane (10 mg/L)
4. Tyramine (25 mg/L)
5. Spermidine (5.0 mg/L)



## Inert Family

“Inertsil” High Performance LC Column

“InertCap” GC Capillary Column


“InertSep” Solid Phase Extraction for pretreatment of analysis



## Examples of sample pretreatment


### Fish samples

### Soy sauce samples




\* Before loading the sample on SPE cartridge, conditioning was carried out with 10 mL of methanol followed by 10 mL of water.

### Extract from fish




### Extract from soy sauce



### Derivatization reaction by DNS-Cl

DNS-Cl reacts with primary and secondary amines. As a result, fluorescent derivatives are generated.



## Contact us

### GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-chome, Shinjuku-ku, Tokyo, 163-1130 Japan

TEL: +81-3 (5323)6620 FAX: +81-3 (5323)6621

### GL Sciences, Inc. USA

4733 Torrance Blvd. Ste 255, Torrance, CA 90503

Tel: (310)265-4424 FAX (310)265-4425

### Distributors Outside of Japan and USA

GL Sciences uses distributors in many countries.

You can find a local distributor in your country in the following url.

<http://www.glsciences.com/products/contact.html>

E-MAIL: [world@glsc.co.jp](mailto:world@glsc.co.jp)

